

# MONOFOCAL IOLs

## PRODUCT PRESENTATION




**ALSIOL**




**ALSIOL VF**



**ALSAFIT**



# 24H Aspherical Correction





Depth of focus

# Day light

Contrast sensitivity

# Twilight

Light scattering control

# Night




## Definition of Sphericity & Asphericity

### SPHERICAL ABERRATION

- When parallel rays pass through the spherical surface, the **central rays** focus on more posteriorly while peripheral rays focus progressively more anterior.
- "Positive" spherical aberration means **peripheral rays** are **bent too much**.
- Specially at dim night conditions, **Spherical aberrations increase** as a result of expanded pupil size


Lens with Spherical Aberration





## Why Aspheric Correction is so crucial?

- Without aspherical correction on IOL, Aberration causes defocusing of improper refracted lights.
- Total spherical aberrations will result
  - Astigmatism-like aberration
  - Light scatterings in mesopic conditions






## Spherical Aberration Changes in life time period

- In healthy human eye, cornea and crystalline lens have their own spheric aberrations which naturally starts and continue to change from young ages to older ages .
- Consequently both cornea and crystalline lens aberrations effects and determines total **spheric aberrations** in life time period.

TOTAL  
Aberrations





Brunette et al, 2003, IOVS

Recent studies of simultaneous measurements of corneal and ocular aberrations in normal eyes showed that there is a natural compensation of the spherical aberration (SA) between the cornea and lens for young subjects. Normal aging tends to disrupt this compensation mechanism because the SA of the lens gradually changes from negative to positive.



# Balance of corneal and internal SA Averagely remains +0.1 μm SA to have depth of focus in between 20s to 40

Balance between corneal and internal SA remains aprx. 0.1μm which stays constant in between ages 20s to 40.\*



\*Adaptive model of the aging emmetropic eye and its changes with accommodation. Rafael Navarro  
Journal of Vision (2014) 14(13):21, 1–17






## Residual SA Comperison

Aspherical IOLs	24H Aspherical Correction							Spherical IOLs
	Tecnis 1	Acrysof	<b>AlsiOL- AlsaFIT</b>	Zeiss CT ASPHINA 509M	Akreos	Softec HD	enVista	Spherical IOL
Aberration correction of lens	-0.27 $\mu\text{m}$	-0.20 $\mu\text{m}$	<b>-0.09 <math>\mu\text{m}</math></b>	0	0	0	0	+0.15
Residual aberration of eye	0	0.07 $\mu\text{m}$	<b>0.18 <math>\mu\text{m}</math></b>	0.27 $\mu\text{m}$	0.27 $\mu\text{m}$	0.27 $\mu\text{m}$	0.27 $\mu\text{m}$	0.42 $\mu\text{m}$

Strong (-) Aspherical IOLs      Mild (-) Aspherical      Neutral Aspherical IOLs



Expectations & SA correction Comparison	Spherical IOLs	Neutral Aspherical IOLs	Well known Strong (-) Aspherical IOLs	ALSIOL 24H Aspherical Correction
Astigmatism like aberration control	-	+	+	+
Light scatterings and blurred vision control in mezopic conditions	-	-	+	+
Contrast sensitivity at dimlight	-	-	-	+
Decentralization tolerance	-	+	-	+
Depth of focus	-	+	-	+
Leaving a weak spherical value about +0.1 µm as the young eye behaves	-	-	-	+
Focus Profile	Defocus	Defocus	NO longitudinal Focus Narrow focus	Longitudinal focus as young eye



## Provides Deeper depth of focus vision at day & night

Leaves weak spherical value about +0.1  $\mu\text{m}$  as the young adaptive eye behaves


Poor depth of focus



Deeper depth of focus



Thanks to 24H asphericity, Alsanza IOLs control:  
Light scatterings, glare and halos at  
**twilight & night**





Light efficiency improves contrast sensitivity at  
**day and night**



# Smart Matrix Material with IPN technology

**MATERIAL**

**SMART  
MATRIX**



## What is Smart Matrix Material?

- Smart Matrix material with IPN technology is unique combination of 2- Hydroxyethylmethacrylate and 2- Oxyethylmethacrylate special monomers both yielding **4 significant features at same time**







# What is IPN technology?

## Interpenetrating Polymer Network

**Smart Matrix with IPN technology** carries 2 different chemical identities simultaneously providing best features of two **IPN (Interpenetrating Polymer Network) technology builds** Ultra Pure monomers with understroyable connections in **3 Dimentional** structure. That creates a strength polymer network protecting IOL from deformation, delamination or dispersion



2-Hydroxyethylmethacrylate  
Hydrophilic




Oxyethylmethacrylate  
Hydrophobic



**Smart Matrix with IPN technology**





## 1- Both Elastic & Strengthened Material

Smart Matrix Material with IPN technology serves 2 different chemical identities with providing different advantages of two


Expansion	
Linear	1.125
Radial	1.125
Telsile g/mm <sup>2</sup>	33
Young's Modulus g/mm <sup>2</sup>	62



ALSIOLs can expand up to 1.125 times more than its original form linearly and radially

Tensile Strength: ALSIOL can resist stress about 33g/mm<sup>2</sup> without any deformation during all procedure

Young's module: Elasticity module



## 1- Both Elastic & Strengthened Material


- **High Elasticity module** provides IOL folding-inserting and most significantly return ability to its original shape and optical performance in the most desirable time frame
- **High tensile Strength** provides IOL **resist against external** stresses while going through injector to capsular bag in order to prevent any damage on its haptics or optic zone

## 2- 3D integrated Hydrophobic Surface

Thanks to IPN technology, IOL's Smart Matrix material has Hydrophobic surface with undestroyable connections


- Increases fibronectin adhesion
- Lowers affinity to Lens Epithelial Cells (LECs)
- Reduces the risk of PCO





### 3- Pure Material 2 times more Pure than others


- **Ultra Pure monomers and Zero polymerization technology;** lowers the residual methacrylic acid monomer content and results in zero ionicity polymer which is 10 to 20X less acid than other 'highest purity' commercial monomers.
- **Zero ionicity of Smart matrix material** provides further degree of quality by eliminating the possibility of calcium phosphate particles appearing over time. And increases **strength** and **biocompatibility** of ALSanza IOLs
- **Supports High optic Performance** In Point Spread Function Test (PSFT). Pure material keeps IOL durable and prevent microwrinklings on optic zone after folding and unfolding total hydrophobic IOLs additionally




### 3- Pure Material

produced in Virgin Glass Molds to higher biocompatibility

- **Virgin glass molds eliminates the contamination** unlike plastic molds that have high risk to release agents, toxic components, lubricants, antioxidants commonly found in polypropylene and polyethylene resins
- These chemical contaminants can be found in competitor materials produced in plastic (polypropylene, polyethylene) molds.






## 4- Smooth IOL Surface

Smart Matrix exhibits **%33 lower surface roughness** than competitors due to its ultra pure monomer polymerization procedures with **zero ionicity** (less acid). That reduces matrix size and surface roughness of material

Designed for;


- Reduced LEC Migration and adhesion to lens optic
- Reduced PCO incidence\*

The amount of surface irregularities is linearly related to the number of inflammatory cells adhering to the optic surface of the IOL and to the rate of LEC migration and PCO incidence\*



\*Analysis of intraocular lens surface properties with atomic force microscopy Marco Lombardo, MD, Maria P. De Santo, PhD, Giuseppe Lombardo, PhD, Riccardo Barberi, PhD, Sebastiano Serrao, MD, PhD

MATERIAL  
SMART  
MATRIX



# FIT HAPTICS

## Arms2FIT

## Fins4FIT





## Expectations from Haptics in Surgery

A- **Responsible for stabilization** of IOLs inside crystalline bag.  
That is a key preventing factor against to rotation and tiltations.

B- Proper unfolding without **any harm on crystalline bag** during procedures

C- **Suspensive role** of haptics is **main protective factor** for desired optical results




D- Preserve **refractive stability with resisting contractive stresses in bag**



## Definifition of Haptic Suspension

- **FIT HAPTICS** are designed to resist tensile and bending stresses which comparely high in upper and lower part of crystalline lens bag
- Crsytalline bag pressures are transmitted from FIT Haptic contact areas to IOL body
- FIT HAPTICS specially designed to behave **suspensive against external stresses in order to preserve stability and desired central optic position.**


When compressed or stretched, a spring gains elastic potential energy.



With suspensive behavior, stability and central optic position is preserved against high tensile and bending stresses




# What is the difference from the Dr's perspective interms of Arms2Fit and Fins4FIT designs?

FIT HAPTICS		
	Arms2Fit	Fins4Fit
	Overall Ø : 13.00 mm Classical model	Overall Ø : 11.00 mm New generation
	2 contact points	4 contact points
	<b>2 Suspensive Arms</b>	<b>4 Suspensive Fins</b>
	<b>Preserved stability and desired central optic position</b>	<b>Preserved stability and desired central optic position</b>



## Arms2Fit



## Arms2Fit

### Preserve stability and desired central optic position

Arms2FIT Haptic design **adds** a new perspective to tradition with its unique design

#### Unique Suspensive additions;

- **‘Narrow junction angle’** improves suspensive movements of Arms to maintain IOL and optical field stability
- **‘Stress Dimmer Curves’** of inner arms **softens transmitted stresses** to support function of Suspensive Narrow junction angle



## Arms2Fit

### Preserve stability and desired central optic position

#### Both Flexy & Strong FIT Arms


- **Minimize** tiltations mainly resulting from contractions, tensile or bending stresses sourced by crystalline bag
- Support **Harmony** between Crystalline Bag & IOL to preserve desired refractive stability

**Overall 13mm distance**




# Arms2FIT


## FIT Arms and Stress dimmer curves prevent reaching of stress to the optical field.



9,5 mm



10 mm



13 mm

## Fins4FIT

4 Suspensive  
corner Fins

2 Wide Aspiration and  
4 Positioning Holes






# Fins4FIT

## Preserve stability and desired central optic position

### Sophisticated Design

- **Suspensive Corner Fins**; behave self-adaptive to fit different sizes of cristalline bag
- **4 Wide OVD Aspiration Holes**; facilitate final aspirations in surgery
- **4 positioning holes**; supports more accurate lens positioning



0° of Fin4FIT haptic angle


# Fins4FIT

## Preserve stability and desired central optic position


### II. Extremely Suspensive

- 4 Suspensive Corner Fins demonstrate suspensive movements against compression stresses. Show high tendency to preserve IOL position and refractive stability with reducing risks of tiltations caused eventual crystalline bag contractions

- 6 Additional holes are also supportive adds to soften transmitted external stresses on to IOL



When compressed or stretched, a spring gains elastic potential energy.



## Fins4FIT

### Preserve stability and desired central optic position



#### III. Harmonically stable and adaptive

- Fins4FIT haptic design enhances ALSAFIT IOLs with self-adaptive feature of corner fins. Corner Fins able to stabilize all IOL either in Small or Large Capsular bag sizes by easy self-adaptation

## Fins4FIT

### Preserve stability and desired central optic position

- Self-adjustable behavior of 4 Suspensive Corner Fins in various Capsular Bag sizes.



- I. Sophisticated Design**
- II. Extremely Suspensive**
- III. Harmonically stable and adaptative**



# DIAMOND CUT EDGE

**DIAMOND CUT  
EDGE**





## Posterior Capsule Opacification (PCO) types

Preventing “after-cataract” reactions is a strenuous fight against nature because lens epithelial cell (LEC) migration and transdifferentiation are nature’s back-up strategies to the loss of lens tissue.

Eyes showing various forms of capsular opacification.



*A, Extensive anterior capsule opacification;*

*B, Mixed type of PCO with fibrous (arrow) and pearl type areas (asterisk);*

*C, Fibrous form of posterior capsule opacification (arrow);*

*D, Proliferative or pearl form of posterior capsule opacification (asterisk);*


*E, Linear posterior capsule opacification*





## PCO still a major hurdle in successful cataract surgery

- PCO results from the growth and **abnormal proliferation of LECs** (Lens Epithelial Cells) on the capsule at the time of cataract surgery.
- These cells **migrate** to the posterior capsule where they approach the central visual axis and cause visual axis obscuration, resulting in dimness of vision.








## 'Edge Effect' Junctional area of Optic & Haptic parts of IOL

- Recent clinical experiences shows that **'edged IOLs reduces the PCO area and severity** compared with an identical **round-edged IOLs'**.
- According to Kerry D. Solomon, MD, OSN Refractive Surgery Board Member, he noted **'edged design slows down or delay the migration,** and that carried over to the design of lens optics."  
"The more square of an edge and the more defined of an edge, the better the delay or the better the barrier"







## Continous Optic Quality

All ALSANZA IOLs have most prominent edges with Diamond Cut Edge design that frames all optical field. More sharpened design of Diamond Cut Edge most efficiently blocks LEC migration than conventional square edge & round edge designs.

- Diamond Cut Edge is developed for Anti-Migration effect for LECs and reduce PCO occurrence



# HIGHEST ABBE NUMBER

Cromatic Aberration  
Control

**CHROMATIC  
ABERRATION**

**HIGHEST  
ABBE  
NUMBER**

## Chromatic Aberration

- **Chromatic Aberration** is a kind of defect of Optic system resulting different wavelengths of **light focus at different distances improperly**
- They are refracted through different angles.
- Chromatic aberration depends on refractive indices of lens. High chromatic aberrations produce a **blurred image with coloured fringes**




**CHROMATIC  
ABERRATION**

**HIGHEST  
ABBE  
NUMBER**

# Chromatic Aberration

Resulting from wavelengths with different refraction index

Longitudinal/Axial  
Chromatic Aberration




- Legend**
- RGB Colour Ray
  - Optical Axis
  - Best Focus Point


## Abbe Number (V - number)

- Chromatic aberration has **inverse ratio** with Abbe Number of IOL

$$V_d = \frac{n_D - 1}{n_F - n_C}$$



IOL material with high chromatic  
aberration  
(low Abbe number)



IOL material with less chromatic  
aberration  
(high Abbe number)

$n_D$ ,  $n_F$  and  $n_C$  are the refractive indices of the material at the wavelengths of the Fraunhofer d- (587.6 nm), F- (486.1nm) and C- (656.3nm) spectral lines. The low dispersion materials have high value of V.

CHROMATIC  
ABERRATION

HIGHEST  
ABBE  
NUMBER

## HIGHEST ABBE NUMBER

- Comparing with others AlsiOLs offers **highest Abbe Number** in market which enables to have strongest Chromatic Aberation Control


	SA correction	Refractive Index	Abbe Number
Alcon Acrysof	-0.20	1.55	37
AMO technis	-0.27	1.47	55
Crystalline Lens	+0.27	1.41	47
<b>ALSIOLs</b>	<b>-0.09</b>	<b>1.46</b>	<b>58</b>

# VIOLET FILTER

## Effective Photoprotection








## Visible Light

The visible spectrum is the portion of the electromagnetic spectrum that is visible to the human eye. A typical human eye will respond to wavelengths from about 390 to 700 nm.






## Violet Filter or Blue blocking ?

Violet Filter reduce retinal exposure to Violet (400-440 nm) light and shorter wavelength including UV radiation.


Blue blocking IOL attenuate blue (440-500) light and shorter wavelengths including UV radiation.





## Why Blue Light Transmission is important

- Blue light has an important role in the **circadian system**, affecting the sleep–wake cycle, mood and cognition.
- Circadian system primarily depends on photopigment **regulation of melanopsin**
- Melanopsin has a peak absorption at **465–480 nm**





## VF: Transmission performance same as Natural lens

- ALSIOL VF is effective for UV and Violet Spectrum photoprotection avoiding total blue blocking
- ALSIOL VF provides %95 of blue light transmission at **480 nm** wavelength that has significant role in the circadian system
- Additionally blue light is **important for the scotopic vision**



## Which chromophore does VF use

- ALSIOL VF has **3- hydroxykyurenine** same photoselective chromophore as natural crystalline lens
- 0,02 % concentration of chromophore do not cause **chromatic disorders** unlike 0,04 % concentration synthetic dye filtered Blue Blocking IOLs does.
- Using natural chromophore and its lower concentrations allow **clearer contrast sensitivity**

**0,02 %**  
3-hydroxykyurenine



**0,04 %**  
Synthetic dye filter chromophore





## VF Filter ...to Arrive Natural Benefit...

- **Mimics Nature**

- Protects retina without blocking essential blue light in order to facilitate circadian system
- Optimal concentration of exact crystalline chromophore provides clearer contrast sensitivity and color perception contrary to 'synthetic yellow filters' does.

- **Offers Sufficient Photoprotection without Deficient Phototransmission**  
For better scotopic vision with blue photoreception

# Manufacturing and Product Quality






## DIOPTER RANGE

ALSIOL, ALSIOL VF, ALSAFIT	FROM -10 TO -0.5 & 32.5 TO 45
ALSIOL, ALSIOL VF, ALSAFIT	FROM -20 TO -10.5 & ABOVE 45
EXTREME PRODUCTION	


Monofocal - Aspheric - FINS4FIT Haptics

**ALSAFIT**



Monofocal - Aspheric - FINS2FIT Haptics

**ALSIOL**



Monofocal - Aspheric - FINS2FIT Haptics - VF

**ALSIOL VF**







## HIGHER DIOPTRIC ACCURACY Well above the Standarts

Alsanza Optic Intellegence standardize **more sensitive dioptric deviation control** than international standards. Manufacturing process ensures to perform **High Dioptric Accuracy** in total range of IOLs with 0.25D very low dioptric deviation value

DIOPTRIC DEVIATION VALUES	
International Standarts: Acceptable max. dioptric deviation value	$\pm 0.4$
<b>ALSANZA Optic Intellegence: Acceptable max. dioptric deviation value</b>	<b><math>\pm 0.25</math></b>





## MTF

- **MTF (Modulation Transfer Funtion)** indicates the ability of an optical system to reproduce (transfer) various levels of details (spatial frequancies) from object to image.
- This is a way to describe the **contrast sensitivity** of an optical system.
- MTF is described as visual preformance of IOLs
- According to international IOLs standards MTF result must be above **0.43 at 100 lpm (lines per milimeter)**



## MTF and Contrast Sensitivity

As line spacing decreases (i.e. **the frequency increases**) on the object, it becomes increasingly difficult for the lens to transfer image contrast. Decreased contrast sensitivity determines low MTF value.





## MTF - ALSANZA Well above the Standards

- **With 0,60 at 100 lpm**, ALSANZA IOLs presents one of highest MTF value in competition
- **Optic quality of ALSANZA IOLs present well above International standards.** Permitted MTF value is established as only 0,43 at 100 lpm for global IOL manufacturers
- ALSANZA executes strict control mechanism to maintain high standards High Standards




Each product's MTF value is individually checked in production line one by one



**24H ASPHERICAL CORRECTION**  
Asphericity for All Day and Night




**MANUFACTURING & OPTIC QUALITY**  
Well above the Standards



**VIOLET FILTER**  
Arrive Natural Benefit



**SMART MATRIX WITH IPN TECHNOLOGY**  
3D structured Privileges



**FIT HAPTICS**  
Preserve stability and desired central optic position



**DIAMOND CUT EDGE**  
Continuous Optic Quality

**HIGHEST ABBE NUMBER**  
Chromatic Aberration Control





**Thank you**