



**ALSAVISC**

Sodium Hyaluronate Viscoelastics

**ALSANZA**  
Medizintechnik und Pharma GmbH

## Visco-Elasticity



Molecular Weight

Surface Tension

DISPERSIVE  
SODIUM HYALURONATE<sup>1</sup>

Diffuse  
3.0%

COHESIVE  
SODIUM HYALURONATE<sup>1</sup>

Cool    Deep  
1.4%    1.8%  
1.6%



- **Molecular Weight: 750.000 Daltons**
- Low surface tension provides surface affinity for endothelium protection
- Higher coating ability for safe procedures


- **Molecular Weight: 3.000.000 Daltons**
- Higher molecular weight & Long HA chains provide easy aspirations
- Shows high elastic behaviour that enables better stress absorption
- Reversible viscosity with shear thinning provides smooth instrument movements and injections

Coating & Protection


Space Maintaining & Shock Stress Absorption

# Visco-Elasticity

**Viscosity of Alsaveics decreases** to facilitate instrument manipulations and smooth injections by ‘thinning’ as shear stress increases.



**Elasticity of Alsaveics** increases with shear stress. Alsaveics have tendency to maintain space and absorb high-compression stresses during phaco with their increasing elastic response.




Elastic response of Alsaveics simultaneously facilitates the return of the viscosity to its original state when shear stress ends.

# Molecular Weight

Cohesive Alsavisc line has **3.000.000Da**

Molecular Weight Sodium Hyaluronate

- Well entanglement with long chains
- Easy aspirations
- Better space maintenance & shock absorption with higher visco-elastic properties



Cool	Deep
1.4%	1.8%
1.6%	


Dispersive Alsavisc Diffuse has **750.000Da**

Molecular Weight Sodium Hyaluronate

- Well dispersion with short chains
- Higher surface affinity to endothelium cells
- Better coating and retention performance in eye during phaco.



Diffuse  
3.0%



Expectations	CCC	Phaco	IOL implantation	Final Asp.
Ease of Injection	<div style="display: flex; gap: 5px;"> <div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div> <div style="background-color: #009688; color: white; padding: 2px; font-size: 8px;">DISPERSIVE SODIUM HYALURONATE</div> </div>			
Creating and maintaining Space	<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>		
Corneal Endothelial Cell & Lens Epithelial Cell Protection - Coating	<div style="background-color: #009688; color: white; padding: 2px; font-size: 8px;">DISPERSIVE SODIUM HYALURONATE</div>	<div style="background-color: #009688; color: white; padding: 2px; font-size: 8px;">DISPERSIVE SODIUM HYALURONATE</div>	<div style="background-color: #009688; color: white; padding: 2px; font-size: 8px;">DISPERSIVE SODIUM HYALURONATE</div>	
Supporting Instrument Manipulation	<div style="display: flex; gap: 5px;"> <div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div> <div style="background-color: #009688; color: white; padding: 2px; font-size: 8px;">DISPERSIVE SODIUM HYALURONATE</div> </div>	<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	
Shock stress absorption		<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	
Deepen and stabilize AC		<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	
Balancing the pressure of posterior and anterior chamber		<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	
Absorbing the pressure created by IOL			<div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div>	
Ease of removal				<div style="display: flex; gap: 5px;"> <div style="background-color: #e91e63; color: white; padding: 2px; font-size: 8px;">COHESIVE SODIUM HYALURONATE</div> <div style="background-color: #009688; color: white; padding: 2px; font-size: 8px;">DISPERSIVE SODIUM HYALURONATE</div> </div>

**COHESIVE**  
SODIUM HYALURONATE

# Cool Deep

		Rheology	Concentration	Molecular Weight	Zero Shear Viscosity	Volume	Origin	Cannula	Osmolality	pH Value	Storage Condition
Cool	<b>Alsavisc 1.4%</b>	Cohesive	14 mg/ml	3 MDa	200.000 mPas	1.1 ml	Biofermented	27 G	300-350 mOsm/kg	6.8-7.6	2-25 °C
	<b>Alsavisc 1.6%</b>	Cohesive	16 mg/ml	3 MDa	400.000 mPas	1.1 ml	Biofermented	27 G	300-350 mOsm/kg	6.8-7.6	2-25 °C
Deep	<b>Alsavisc 1.8%</b>	Cohesive	18 mg/ml	3 MDa	600.000 mPas	1.1 ml	Biofermented	27 G	300-350 mOsm/kg	6.8-7.6	2-25 °C

Table taken from reference 3

**DISPERSIVE**  
SODIUM HYALURONATE

# Diffuse

		Rheology	Concentration	Molecular Weight	Zero Shear Viscosity	Volume	Origin	Cannula	Osmolality	pH Value	Storage Condition
Diffuse	<b>Alsavisc 3.0%</b>	Dispersive	30 mg/ml	0.75 MDa	30.000 mPas	1.1 ml	Biofermented	25 G	300-350 mOsm/kg	6.8-7.6	2-25 °C

Table taken from reference 3

References:

- 1.Review of Ophthalmology Understanding and Using the Full Spectrum of OVDs Devgan, MD, FACS, Sun Valley, Calif. PUBLISHED 2009
- 2.Curr Opin Ophthalmol. 2008 Jan;19(1):50-4
- 3.Data on file.

